

MANUAL

MQ Creator Edition

Musiquence Technologies, Lda


*Empowering creators with AI-driven, science-backed, no-code tools for
interactive, user-centered solutions.*

Est Conde Carvalhal 177a, 9060-014, São Gonçalo Funchal,
Madeira Island - Funchal

NIF: 517244632

Email: info@musiquencetechnologies.pt

Tel: +351962155521



MUSIQUENCE TECHNOLOGIES
be creative!

Index

Index.....	1
Musiquence Technologies, LDA	2
Who are we?.....	2
1. Introduction.....	3
1.1. Purpose of the manual	3
1.2. Scope of application.....	3
2. First steps.....	4
2.1. System requirements.....	4
2.2. Installation instructions	4
2.3. Troubleshooting:	4
3. User interface overview.....	5
3.1. Home page.....	5
4. How to use MQ Creator Edition tools.....	6
4.1. Activities	6
4.2. Editor	7
4.3. Game	9
5. How to use the MQ AI Toolset integrated in MQ Creator?	10
6. Predefined shortcuts in MQ-Editor	13

Musiquence Technologies, LDA

Who are we?

Musiquence Technologies, LDA is the first start-up based in the Autonomous Region of Madeira, created to transfer academic knowledge from the world of research to professional practice.

Our mission is to transform the next generation of artificial intelligence (AI) solutions applied to healthcare and education, combining scientific innovation with a genuine focus on people's well-being.

We seek to improve quality of life by connecting clinical and educational services with their communities, promoting knowledge and care in an integrated and accessible way.

Our products are aimed at teachers, therapists, and caregivers, enabling the customization of educational and therapeutic activities. By automating and simplifying the creative process, we enable quick and effective adaptations without the need for advanced technical knowledge.

Our technology stands out for its intuitive, easy-to-use use and is developed to meet the individual needs of students and users in a practical and efficient way.

Our solutions empower educators and healthcare professionals to implement cognitive stimulation activities based on the use of elements such as music and reminiscence, using the Musiquence Engine ©.

Our proprietary, patent-pending technology, the result of years of research and the award-winning scientific and innovation awards, implements gamified activities with audio and interactive assistance.

1. Introduction

1.1. Purpose of the manual

This manual provides comprehensive instructions on how to use Musiquence Creator Edition—a standalone application where you can create, save, edit, and play your own activities. You can create activities for cognitive stimulation, education, exergames, and more!

This manual includes instructions for installing, configuring and using all integrated tools, ensuring an optimized workflow for healthcare professionals and Musiquence users.

1.2. Scope of application

This manual focuses on Musiquence Creator Edition and does not cover third-party applications or external hardware components.

Target audience

This guide is intended to:

- New users looking to understand the toolset's functionality.
- Existing users looking for advanced guidance on automation and personalization.
- Healthcare professionals who use Musiquence for therapeutic activities.
- Educators who use Musiquence for educational activities.

Basic computer and software proficiency is assumed.

2. First steps

2.1. System requirements

Before proceeding with the installation, verify that your system meets the following requirements:

- Operating System: Windows 10/11
- Minimum RAM: 8 GB
- Minimum disk space: 1 GB
- Internet connection required for online AI models

2.2. Installation instructions

Follow these steps to install Musiquence Creator Edition:

1. Download and install ORI (OpenRehab Initiative) from <https://openrehab.arditi.pt/download-app>.
2. If you don't have an account, sign up for one. Otherwise (or after registering), log in with your credentials.
3. Refresh the apps list by clicking the circular arrows icon (↻) next to the profile icon.
4. Find "MQ Creator Edition" in the list of available applications.
5. Click the three dots (:), select the latest version, and install the toolchain.
6. After installation, click the Start button (▶). If prompted, select the latest version to open the Musiquence Creator Edition app.

2.3. Troubleshooting:

Question	Solution
Unable to log in	Make sure you are using the correct credentials in ORI. Additionally, you must agree to the terms and conditions before registering or logging in. If the problem persists, please speak to ORI administrators.
Cannot find "MQ Creator Edition"	Make sure to refresh the list by clicking the circular arrows icon (↻).

If the problem persists, speak to the administrators of contactinfo@musiquencetechnologies.com.

3. User interface overview

3.1. Home page

The home page provides access to all tools:

- *Language*:select the flag that corresponds to the language in which the application is used.
- *Definitions*:select the technology(ies) and means of interaction to be used.
- *Editor*:create activities or edit pre-made activities.
- *User*:fill in if you want to add a user to store data.
- *To load*:load a pre-made activity to play it.
- *To go out*:close the application.




Figure 1. MQ Creator Edition home page.

4. How to use MQ Creator Edition tools

4.1. Activities

The Musiquence platform offers a range of customizable activities designed to foster user creativity. These include:

- **Quiz**—according to the instruction, the player must select the correct answer(s).
- **Advanced Quiz**—variation of the quiz activity, with greater structuring flexibility.
- **Search**—exploratory research activity, in which the player must find the target stimuli, using a virtual magnifying glass.
- **Association**—association activity where the player must categorize or associate the stimuli, according to the instructions presented.
- **Break**—break activity, ideal for interspersing between other activities, giving the player time to prepare.

In addition to the ease of customizing the content, these activities can be played on different technologies, such as computers, tablets, motion sensors, interactive tables, augmented reality and virtual reality through **MQ player**. The high level of customization allows the user to intervene **personalized and differentiated**, enhancing its impact.




Figure 2. Example of using activities created with MQ Creator Edition in MQ player, using a virtual reality headset.

4.2. Editor

The Editor page brings together all the tools needed to create new activities or edit pre-made ones. Either way, this is a flexible tool that allows users to unleash their creativity. All stimuli—instructions, responses, visual and auditory stimuli—are customizable.

In the top left corner, you'll find several buttons to load a pre-made activity, create an activity from scratch, or save your changes.




Figure 3. Editor page – buttons to load, start or save an activity, or autosave (5 to 5 minutes)

In the "settings" button, you can select the game's general settings—that is, settings that will apply to all activities created in the same file. For example, you can set the font type and size, a general background image and/or sound (which will apply to all slides), display a progress bar during the game, or even the presence or absence of metrics.




Figure 4. “settings” button, on the Editor page.

In the upper right corner, you'll find a button for each of the available activities—Quiz, Quiz Advanced, Search, Association, and Break. Whenever you click one of these buttons, the selected activity will appear in the left bar, in the order selected. However, you can change this order using the arrows next to the slide, or delete.




Figure 5. Editor page features.

All content is customizable, flexible, and easily adaptable. You can preview the edited activity by clicking "preview"—in this mode, you can drag and reposition the included elements, as well as edit the answer options.

Did you know?

Musiquence has a unique feature: it uses musical feedback to guide behavior toward the correct response.

If you want to use this differentiating feature, simply:

- 1) Add background music to your activity.
- 2) Click on "settings".
- 3) Select the distortion you want to use: pitch+ rhythm, pitch or rhythm.

Reminder:

- Musiquence needs to distinguish the correct answer options. Don't forget to select them by clicking the correct (✓) button.
- **Always keep** your activities before proceeding to game mode or closing the application.

4.3. Game

To start game mode, you must have already created activities. You can access game mode in three ways:

- 1) On the home screen, click "Upload" and select the desired activity.
- 2) On the home screen, fill in the user field and select the desired activity
- 3) On the Editor page, click "Game" and select the desired activity.

Then, just play!

Notice:

If you want to save performance data and associate it with a profile, don't forget to fill in the user information (Figure 1 – User). You should always use the same name (User ID) and can use the "Experimental Condition" field to identify the session, and enable "Activate Logging of Performance Metrics (.csv)" (Figure 4).

5. How to use the MQ AI Toolset integrated in MQ Creator?

The Musiquence platform allows you to create activities, images, and even text-to-speech features (patent-pending) using the MQ AI Toolset. However, some of these features can be used directly from the MQ Creator editor. To do this, you need the Toolset to be running in the background and use the following features:

1. If you want to add an AI-generated image, simply click on the “Robot” button that is associated with the images (Figure 6).




Figure 6. “Robot” buttons to create images using AI.

2. If you want to use voice instructions using AI-synthesized voice, just click the button next to the instructions field (Figure 7).




Figure 7. "Robot" buttons to activate AI-synthesized voice.

3. Finally, if you want to create multiple activities instantly, just click on the "Robot" button in the menu (Figure 8).




Figure 8. "Robot" buttons to create activities using AI.

It's important to note that the quality of outputs depends on the type of AI model being used at the time of creation. However, users are free to use the MQ Editor to adapt or modify activities according to their preferences.

Notes:

If you want to use MQ-AI, simply connect the MQ Toolset to the ORI. The MQ platform takes care of the rest. 😊

Whenever you click a Robot button, a text field will appear. Just type whatever you want in there. For example, if you want sports activities, you could enter the following text:

*"Create 6 QUIZ activities on generic facts
of Portuguese football. It also introduces images"*


6. Predefined shortcuts in MQ-Editor

MQ-Editor has a set of predefined shortcuts that help you speed up the activity creation process. The following table lists the shortcuts available in the editor.

Shortcut	Functionality
Generic Actions	
CTRL + Z	Undo
CTRL + SHIFT + Y	Redo
CTRL + S	Fast Recording
CTRL + V	<i>Copy/Paste</i>
CTRL + N	New Project
CTRL + L	Load activity
CTRL + ↑	Navigate to previous slide
CTRL + ↓	Navigate to the next slide
CTRL+I	Create activities (with AI)
Activity Creation	
CTRL+1	Quiz Activity
CTRL+2	Advanced Quiz Activity
CTRL+3	Search Activity
CTRL+ 4	Association Activity
CTRL+5	Coffee Break Activity

Prepared by: Monica Spínola and Luís Ferreira

Revision Date: July 9, 2025



MUSIQUENCE TECHNOLOGIES
be creative!